

- Edits to the noise marking GUI
 - Adding peak marking and heartrate
 - The heartrate is now calculated by a minimum of 10s even if the window currently displayed is shorter than 10s
 - Add hotkeys (1-9) to the noise marking
 - Fixed a bug in the scaling that caused it not to reset when checkbox unchecked
 - Made the annealing, peak finding, and templating run in the background to a file after it has been marked
- Edited the Templating GUI
 - Added in a shift to the PPG so it is accurately displayed as to the right of the ECG
 - This doesn't work if the PPG signal is shifted to the right due to recording error(?)
 - Added right and left to the end of the displayed PPG signal so you can better mark the onset and end
 - Added 9 markers that the user can place around the PPG/ECG
- CHAOS data combination
 - I made an error where I did not combine CHAOS datasets that had a gap - this was fixed and now even if there is a gap between two CHAOS files of the same Patient ID and Room number, it is combined up until 8 hours
 - Got 248 chaos files from the server and ran a script to convert them to bin files